


CONFECTIETABEL PERFECT-FIT

VenoTrain® micro | VenoTrain® business
VenoTrain® soft | VenoTrain® soft S

BAUERFEIND.COM




Omvang	XS*	S		M		L		XL			
	Normaal	Normaal	Plus	Normaal	Plus	Normaal	Plus	Normaal	Plus		
cG siliconenboord wijd*		51–57	59–66	56–63	62–69	60–67	66–74	63–71	71–80		
cG siliconenboord normaal	44–50	45–51	52–59	49–56	55–62	53–60	58–66	55–63	62–71		
cG		45–57	52–66	49–63	55–69	53–67	58–74	55–71	62–80		
cF	39–50	40–52	43–55	43–55	46–60	45–59	50–64	49–63	52–68		
cE	32–38	33–40	36–44	35–43	39–49	38–48	42–53	42–53	47–58		
cD	28–31	29–35	32–38	32–38	35–42	34–42	37–45	37–45	41–49		
cC	29–33	31–37	34–40	34–40	37–44	37–45	39–47	40–48	43–51		
cB1	23–26	25–29	27–31	27–32	29–34	29–35	31–37	32–38	34–40		
cB	18–20	20–22		22–25		25–28		28–31			
cY	27–31	29–33		31–35		34–40		37–44			
cA	17–22	19–24		22–27		25–31		28–35			
Lengte	short		long		Type	*VenoTrain micro AD, XS alleen verkrijgbaar in Ccl. 2, normaal, caramel, gesloten teenstuk VenoTrain soft AD, XS alleen verkrijgbaar in Ccl. 2, normaal, caramel, gesloten teenstuk (zonder siliconenboord)					
lG	61–70		70–80		AG, AT						
lD	35–40		40–46		AD						

Gegevens in het maatsysteem in cm: l staat voor lengte · c staat voor omvang

CONFECTIETABEL

VenoTrain® delight

BAUERFEIND.COM


Omvang	I	II	III	IV	V	VI
cD	28 - 34	30-37	32-39	34-41	36-43	38-46
cC	29-35	31-38	33-40	35-42	37-44	39-47
cB1	23-28	25-30	27-32	29-35	31-37	33-39
cB	19-20,5	20,5-22	22-24	24-26	26-28,5	28,5-31
cY	25-29	26-31	28-33	30-35	32-38	35-41
cA	19-22	20-23	22-25	24-27	26-30	28-33

Kijk voor maatwerkuitvoeringen in de actuele productcatalogus.

Lengte	short	long	Type
ℓD	35-39	40-44	AD