

# mediven elegance®

round-knitted

<b>Product</b>	<b>mediven elegance</b> All you see is beauty
<b>Material</b>	light
<b>CCL</b>	1, 2
<b>Lengths</b>	Standard: AD, AG, AT, AT / U
<b>Standard colours</b>	Black, Anthracite, Navy, Cashmere, Caramel, Rosé, Bronze, Beige, medi Magenta*, White*
<b>Sizes</b>	I – VII
<b>Standard /made-to-measure</b>	Standard and made-to-measure
<b>Foot part</b>	Open and closed toe, universal foot size, shoe sizes 36 – 43
<b>Topband / attachment</b>	Silicone dot topband (NoB), Silicone dot topband motif (NoM), Plain silicone topband rose (FBR), Silicone micro dot topband sensitive (MBs). All topbands are standard styles.
<b>Special versions</b>	Possible
<b>Key words</b>	Transparent, Soft Elastic foot tip, wide and soft cuff, smaller heel for wrinkle-free fit.
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\* Delivery period up to 10 working days.

# Measuring chart

mediven elegance®

Standard size seamless		I	II	III	IV	V	VI	VII	
AG with topband normal	cG	43-48	45-52	49-56	53-60	56-64	60-68	64-72	
AG with topband wide	cG	49-57	53-62	57-67	61-72	65-77	69-81	73-85	
<p>l=Lengths c=Circumferences</p>	cG	43-57	45-62	49-67	53-72	56-77	60-81	64-85	
	cF	39-52	41-56	44-60	47-65	50-69	53-73	57-78	
	cE	30-37	33-40	35-43	37-45	39-48	41-51	42-53	
	cD	27-33	29-36	32-39	34-42	36-45	38-48	40-50	
	cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51	
	cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41	
	cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32	
	cY	26-31	28-33	29-35	31-37	32-38	33-40	34-42	
	cA	17-22	19-24	21-26	23-29	25-32	27-34	28-35	
	<p>Lengths</p>	lG normal	72-83 cm		<p>Foot</p>	Highly elastic universal foot size			
		lG petite	62-71 cm						
		lD normal	39-44 cm						
lD petite		34-38 cm							

# mediven® for men

round-knitted


<b>Product</b>	<b>mediven for men</b> The elegant compression sock for men
<b>Material</b>	light
<b>CCL</b>	1,2
<b>Lengths</b>	AD
<b>Colours</b>	Black, Anthracite, Navy, Brown
<b>Sizes</b>	I – VII
<b>Standard / made-to-measure</b>	Standard
<b>Foot part</b>	Closed toe, shoe sizes 39 – 46
<b>Key words</b>	Anatomical toe design, right / left markings, Business Rib Design, Clima Fresh, Added Lanolin*.
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\* Natural substance used for treating dry and chapped skin, has a broad range of applications, for example in child care and baby care. The caring protective effect remains even after washing.

# Measuring chart

mediven® for men

**Standard size seamless**      I      II      III      IV      V      VI      VII



<b>Circumferences in cm</b>	<b>cD</b>	27-33	29-36	32-39	34-42	36-45	38-48	40-50
	<b>cC</b>	28-34	30-37	33-40	35-43	37-46	39-49	41-51
	<b>cB¹</b>	23-27	24-29	26-32	29-35	31-37	33-39	35-41
	<b>cB</b>	18-20	20-22	22-24	24-26	26-28	28-30	30-32
	<b>cY</b>	26-31	28-33	29-35	31-37	32-38	33-40	34-42
	<b>cA</b>	17-22	19-24	21-26	23-29	25-32	27-34	28-35

<b>Lengths</b>	<b>mediven for men</b>	
	<b>ℓD normal</b>	39-44 cm
	<b>ℓD petite</b>	34-38 cm

# mediven® comfort

round-knitted


<b>Product</b>	<b>mediven comfort</b> All you feel is softness
<b>Material</b>	medium
<b>CCL</b>	1, 2
<b>Lengths</b>	Standard: AD, AG, AT
<b>Standard colours</b>	Black, Anthracite, Navy, Cashmere, Caramel, Rosé*, Bronze*, Beige*
<b>Sizes</b>	I – VII
<b>Standard / made-to-measure</b>	Standard and made-to-measure
<b>Foot part</b>	Open and closed toe, universal foot size, shoe sizes 36 – 43, double-knitted comfort sole
<b>Topband / attachment</b>	Silicone dot topband (NoB), Silicone dot topband motif (NoM), Plain silicone topband rose (FBR), Silicone micro dot topband sensitive (MBs). All topbands are standard styles.
<b>Special versions</b>	Possible
<b>Key words</b>	Soft material, double-knitted comfort sole, added Lanolin**.
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\* Delivery period up to 10 working days.

\*\* Natural substance used for treating dry and chapped skin, has a broad range of applications, for example in child care and baby care.  
The caring protective effect remains even after washing.

# Measuring chart

mediven® comfort

Standard size seamless		I	II	III	IV	V	VI	VII	
AG with topband normal	cG	43-48	45-52	49-56	53-60	56-64	60-68	64-72	
AG with topband wide	cG	49-57	53-62	57-67	61-72	65-77	69-81	73-85	
 <p>ℓ=Lengths    c=Circumferences</p>	cG	43-57	45-62	49-67	53-72	56-77	60-81	64-85	
	cF	39-52	41-56	44-60	47-65	50-69	53-73	57-78	
	cE	30-37	33-40	35-43	37-45	39-48	41-51	42-53	
	cD	27-33	29-36	32-39	34-42	36-45	38-48	40-50	
	cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51	
	cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41	
	cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32	
	cY	26-31	28-33	29-35	31-37	32-38	33-40	34-42	
	cA	17-22	19-24	21-26	23-29	25-32	27-34	28-35	
	ℓB <sup>1</sup>	cB <sup>1</sup>	<p><b>Foot</b> Highly elastic universal foot size Double-knitted comfort sole</p>						
	ℓG normal	72-83 cm							
	ℓG petite	62-71 cm							
	ℓD normal	39-44 cm							
ℓB	cB	<p><b>Lengths</b></p>	ℓD petite	34-38 cm					
ℓC	cC								
ℓD	cD								
ℓE	cE								
ℓF	cF								
ℓG	cG								

# mediven plus®

round-knitted

	<b>mediven plus</b> All you experience is reliability
<b>Material</b>	medium
<b>CCL</b>	1, 2, 3
<b>Lengths</b>	Standard: AD, AG, AT, AT / U, AT / H, AT / E
<b>Standard colours</b>	Black, Anthracite, Navy, Cashmere, Caramel, Rosé, Bronze, Beige, medi Magenta*, White*
<b>Sizes</b>	I – VII
<b>Standard / made-to-measure</b>	Standard and made-to-measure
<b>Foot part</b>	Open and closed toe, universal foot size, shoe sizes 36 – 43
<b>Topband / attachment</b>	Silicone dot topband (NoB), Silicone dot topband motif (NoM), Plain silicone topband rose (FBR), Silicone micro dot topband sensitive (MBs). All topbands are standard styles.
<b>Special versions</b>	Possible
<b>Key words</b>	Hard-wearing and opaque, plus heel, optimised toe, ribbed sole reinforcement, 5 cm wide cuff.
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\* Delivery period up to 10 working days.

# Measuring chart

mediven plus®

Standard size seamless		I	II	III	IV	V	VI	VII	
AG with topband normal	cG	43-48	45-52	49-56	53-60	56-64	60-68	64-72	
AG with topband wide	cG	49-57	53-62	57-67	61-72	65-77	69-81	73-85	
<p>l=Lengths    c=Circumferences</p>	cG	43-57	45-62	49-67	53-72	56-77	60-81	64-85	
	cF	39-52	41-56	44-60	47-65	50-69	53-73	57-78	
	cE	30-37	33-40	35-43	37-45	39-48	41-51	42-53	
	cD	27-33	29-36	32-39	34-42	36-45	38-48	40-50	
	cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51	
	cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41	
	cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32	
	cY	26-31	28-33	29-35	31-37	32-38	33-40	34-42	
	cA	17-22	19-24	21-26	23-29	25-32	27-34	28-35	
	<b>Lengths</b>	lG normal	72-83 cm						
lG petite		62-71 cm							
lD normal		39-44 cm							
lD petite		34-38 cm							
<b>Maxi-panty top</b>	<b>AT mediven plus</b> We recommend a maxi-panty top above an abdominal circumference of								
	Size I-III	T = from 95 cm			H = from 110 cm				
	Size IV-V	T = from 120 cm			H = from 135 cm				
	Size VI-VII	T = from 140 cm			H = from 145 cm				
<b>Foot</b>	Highly elastic universal foot size for mediven plus								

mediven active®


round-knitted

<b>Product</b>	<b>mediven active</b> The hard wearing compression sock for men
<b>Material</b>	strong
<b>CCL</b>	1,2
<b>Lengths</b>	AD
<b>Colours</b>	Black, Anthracite, Navy, Brown, White
<b>Sizes</b>	I – VII
<b>Standard / made-to-measure</b>	Standard
<b>Foot part</b>	Closed toe, shoe sizes 38 – 48, comfort sole, seam you can't feel at the toe
<b>Key words</b>	Fine-ribbed design, comfort sole with reinforcement, soft and wide cuff, very robust knitted fabric.
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

# Measuring chart

mediven active®

## Standard size seamless



	I	II	III	IV	V	VI	VII
<b>cD</b>	27-33	29-36	32-39	34-42	36-45	38-48	40-50
<b>cC</b>	28-34	30-37	33-40	35-43	37-46	39-49	41-51
<b>cB<sup>1</sup></b>	23-27	24-29	26-32	29-35	31-37	33-39	35-41
<b>cB</b>	18-20	20-22	22-24	24-26	26-28	28-30	30-32
<b>cY</b>	26-31	28-33	29-35	31-37	32-38	33-40	34-42
<b>cA</b>	17-22	19-24	21-26	23-29	25-32	27-34	28-35

Lengths	mediven active	
ℓD normal	43-49	cm
ℓD petite	37-42	cm

# mediven® forte


round-knitted

	<b>mediven forte</b> All you feel is protection
<b>Material</b>	strong
<b>CCL</b>	2, 3
<b>Lengths</b>	AD, AG, AT
<b>Standard colours</b>	Black, Anthracite, Navy, Cashmere, Caramel, Rosé*, Bronze*, Beige*
<b>Sizes</b>	I – VII
<b>Standard / made-to-measure</b>	Standard and made-to-measure
<b>Foot part</b>	Open and closed toe      Small foot up to 23 cm $\hat{=}$ shoe sizes 36 – 38 Large foot from 23 cm $\hat{=}$ shoe sizes 38 – 43
<b>Topband / attachment</b>	Silicone dot topband (NoB), Silicone dot topband motif (NoM), Plain silicone topband rose (FBR), Silicone micro dot topband sensitive (MBs). All topbands are standard styles.
<b>Special versions</b>	Possible
<b>Key words</b>	Extremely hard-wearing, seamless material, opaque, Clima Fresh over the whole leg, stable soft cuff at the AD border.
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\*Delivery period up to 10 working days.

# Measuring chart

mediven® forte

Standard size seamless		I	II	III	IV	V	VI	VII
AG with topband normal	cG	44-47	46-50	50-54	53-58	56-61	60-65	64-70
AG with topband wide	cG	48-53	51-57	55-61	59-65	62-70	66-74	71-80
 <p>l=Lengths c=Circumferences</p>	cG	44-53	46-57	50-61	53-65	56-70	60-74	64-80
	cF	40-47	42-51	45-55	48-58	51-62	54-67	58-72
	cE	30-37	33-40	36-42	38-44	40-46	42-48	44-50
	cD	28-33	30-35	32-37	34-39	36-41	38-43	40-45
	cC	29-34	31-36	33-38	35-40	37-42	39-44	41-46
	cB <sup>1</sup>	24-26	26-28	28-31	30-33	32-35	34-37	36-39
	cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
	cY	28-30	29-32	31-34	33-36	34-38	36-40	37-41
	cA	18-21	20-23	22-25	24-27	26-29	28-31	30-33
	Lengths	lG normal	72-83 cm					
lG petite		62-71 cm						
lD normal		38-43 cm						
lD petite		33-37 cm						
Foot	Small foot up to 23 cm ≙ shoe sizes 36-38							
	Large foot from 23 cm ≙ shoe sizes 38-43							

# mediven® angio


round-knitted

<b>Produkt</b>	<b>mediven angio</b> Safety in treatment – down to the toes.
<b>Material</b>	medium
<b>CCL</b>	1, 2
<b>Lengths (Standard &amp; MTM*)</b>	AD
<b>Standard colours</b>	Black, Caramel
<b>Sizes</b>	I – VII
<b>Foot part</b>	Closed toe: small foot up to 25 cm $\hat{=}$ shoe size 36 – 40 large foot from 25 cm $\hat{=}$ shoe size 40 – 46
<b>Topband</b>	Tricot border
<b>Special versions</b>	Not possible
<b>Key words</b>	Suitable for concomitant mild to moderate peripheral artery disease (PAD) and / or diabetes mellitus, thanks to Merino-Plush-Padding in the entire foot area up to above the ankles, extra long and flat toe seam, solid ribbed structure, 90°-heel, extra wide cuff, safety proven in clinical trial
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide, Merino wool

\* Delivery period up to 10 working days.

Standard sizes	I	II	III	IV	V	VI	VII
----------------	---	----	-----	----	---	----	-----

ℓ = lengths c = circumferences



Circumferences in cm	cD	27 – 33	29 – 36	32 – 39	34 – 42	36 – 45	38 – 48	40 – 50
	cC	28 – 34	30 – 37	33 – 40	35 – 43	37 – 46	39 – 49	41 – 51
	cB1	23 – 27	24 – 29	26 – 32	29 – 35	31 – 37	33 – 39	35 – 41
	cB	18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32
	cY	26 – 31	28 – 33	29 – 35	31 – 37	32 – 38	33 – 40	34 – 42
	cA	17 – 22	19 – 24	21 – 26	23 – 29	25 – 32	27 – 34	28 – 35

Lengths	Product lengths	
	ℓD regular	39 – 44 cm
ℓD petite	34 – 38 cm	

Foot lengths	
small foot	closed toe, up to 25 cm ≙ shoe sizes 36 – 40
large foot	closed toe, from 25 cm ≙ shoe sizes 40 – 46

# mediven harmony®

round-knitted

<b>Product</b>	<b>mediven harmony</b> Round-knitted armsleeve, flat-knitted hand part
<b>CCL</b>	1, 2
<b>Lengths</b>	Standard length: AG, CG (regular and extra wide) Short length: AG, CG with topband (regular and extra wide) AC <sup>1</sup> hand piece or glove
<b>Colours</b>	Black*, Caramel, Sand
<b>Sizes</b>	I – VIII, armsleeve II – VII, hand part / fingers
<b>Standard / made-to-measure</b>	Standard
<b>Topband / attachment</b>	Narrow napped silicone topband, shouldercap
<b>Key words</b>	Comfortable and effective, double ellipses on the point and in the crease of the elbow
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\*Black: minimum order quantity of 20 pieces per size and variant.

Armsleeve	Size I		Size II		Size III	
Circumference (cm)	Regular	Extra wide	Regular	Extra wide	Regular	Extra wide
cG	25,5 – 28,0	28,5 – 31,0	28,0 – 30,5	30,5 – 33,5	30,0 – 33,0	33,0 – 36,0
cF	23,5 – 25,5	25,5 – 28,0	25,0 – 27,5	27,5 – 30,0	27,5 – 30,0	30,0 – 33,0
cE	23,0 – 25,0	24,0 – 26,0	24,5 – 26,5	26,5 – 29,0	26,0 – 28,5	28,5 – 31,0
cD	22,0 – 24,0	24,0 – 26,0	23,0 – 25,0	25,0 – 27,5	24,0 – 26,0	26,0 – 28,5
cC 1 = 100% **	15,0 – 16,5	15,0 – 16,5	16,5 – 18,0	16,5 – 18,0	18,0 – 19,5	18,0 – 19,5
cC	13,5 – 14,5	13,5 – 14,5	14,5 – 15,5	14,5 – 15,5	16,0 – 17,5	16,0 – 17,5
cA	15,0 – 18,0	15,0 – 18,0	16,0 – 19,0	16,0 – 19,0	17,0 – 20,0	17,0 – 20,0
Armsleeve	Size IV					
Circumference (cm)	Regular	Extra wide				
cG	32,0 – 35,0	35,0 – 38,5				
cF	29,5 – 32,0	32,0 – 35,0				
cE	27,5 – 30,0	30,0 – 33,0				
cD	25,0 – 27,5	28,0 – 30,5				
cC 1 = 100% **	19,5 – 21,0	19,5 – 21,0				
cC	17,5 – 19,0	17,5 – 19,0				
cA	19,0 – 22,0	19,0 – 22,0				

Armsleeve	Size V		Size VI		Size VII	
Circumference (cm)	Regular	Extra wide	Regular	Extra wide	Regular	Extra wide
cG	33,5 – 36,5	37,0 – 40,5	36,0 – 39,5	39,5 – 43,0	38,0 – 41,5	41,5 – 45,5
cF	31,0 – 34,0	34,0 – 37,0	33,0 – 36,0	36,0 – 39,5	35,0 – 38,5	38,5 – 42,0
cE	29,0 – 31,5	31,5 – 34,5	30,5 – 33,5	33,5 – 36,5	32,5 – 35,5	35,5 – 39,0
cD	27,0 – 29,5	29,5 – 32,0	28,5 – 31,0	31,0 – 34,0	30,5 – 33,5	33,5 – 36,5
cC 1 = 100% **	21,0 – 23,0	21,0 – 23,0	23,0 – 25,0	23,0 – 25,0	25,0 – 27,5	25,0 – 27,5
cC	19,0 – 20,5	19,0 – 20,5	20,5 – 22,5	20,5 – 22,5	22,5 – 24,5	22,5 – 24,5
cA	20,0 – 23,0	20,0 – 23,0	22,0 – 25,0	22,0 – 25,0	24,0 – 27,0	24,0 – 27,0

Armsleeve	Size VIII***	
Circumference (cm)	Regular	Extra wide
cG	40,5 – 44,5	45,0 – 49,5
cF	37,5 – 41,0	41,5 – 45,5
cE	35,0 – 38,5	38,5 – 42,0
cD	33,0 – 36,0	36,0 – 39,5
cC 1 = 100% **	27,5 – 30,0	27,5 – 30,0
cC	25,0 – 27,5	25,0 – 27,5
cA	26,0 – 29,0	26,0 – 29,0

Length (cm)	Standard length	Short length
	41 – 47	36 – 40

\*\* Measurement CC1 is 5 cm above CC

\*\*\* Size VIII: Minimum order quantity of 20 pieces per colour and variant; delivery period approx. 4 weeks

## Measuring chart

hand part

### Variants

Seamless AC1 hand piece  
without fingers, with thumb base

Seamless AC1 glove  
with compressive fingers



Length / Hand part  
(cm),  
all sizes

AB	3,5 – 4,5
AC	10,0 – 12,0
AC1	15,0 – 16,0

### Circumference / Hand part (cm)

Sizes	II	III	IV	V	VI	VII
CC 1 = 100 % *	15,0 – 18,0	17,0 – 20,0	19,0 – 21,5	21,0 – 23,5	23,0 – 26,0	25,5 – 30,0
cC	13,5 – 15,5	15,0 – 18,0	17,0 – 20,0	19,0 – 21,0	20,0 – 23,5	22,5 – 27,5
cB	15,5 – 20,0	17,5 – 21,0	19,5 – 23,0	20,5 – 24,0	22,5 – 26,0	24,5 – 30,0
cA	15,0 – 19,0	17,0 – 20,0	19,0 – 22,0	20,0 – 23,0	22,0 – 25,0	24,0 – 29,0



### Length / Fingers (cm)

Thumb	3,5 – 4,5
Index finger	4,5 – 5,5
Middle finger	5,5 – 6,5
Ring finger	4,5 – 5,5
Little finger	3,5 – 4,5


\* Measurement CC1 is 5 cm above CC

## Circumference / Fingers (cm) / Sizes

Sizes		II	III	IV	V	VI	VII
Thumb	Bottom cX	7,0	7,5	8,0	8,5	9,0	10,0
	Top cZ	6,0	6,5	6,5	7,0	7,5	8,0
Index finger	Bottom cX	6,5	7,0	7,5	8,0	8,5	9,5
	Top cZ	5,0	5,5	6,0	6,0	6,5	7,0
Middle finger	Bottom cX	6,5	7,0	7,5	8,0	8,5	9,5
	Top cZ	5,0	5,5	6,0	6,0	6,5	7,0
Ring finger	Bottom cX	6,0	6,5	7,0	7,0	7,5	9,0
	Top cZ	5,0	5,0	5,5	5,5	6,0	6,5
Little finger	Bottom cX	5,5	6,0	6,0	6,5	7,0	8,0
	Top cZ	4,5	4,5	5,0	5,0	5,5	6,0



Circumference CZ is measured at the nailbed.  
Circumference CX is measured at the base of the finger.



### Measurement mediven harmony

\* Length CG is measured on the inside of the arm from the wrist to the anterior axillary fold.

# mediven® 550 leg

flat-knitted

<b>Product</b>	<b>mediven 550 leg*</b> Strong. Very high wall stability. Reliable hold.
<b>CCL</b>	2
<b>Lengths</b>	2 lengths for AD / calf and AD / calf petite
<b>Standard colour</b>	Caramel
<b>Sizes</b>	IV – XIV
<b>Standard / made-to-measure</b>	Standard and made-to-measure
<b>Foot part</b>	Open toe straight ending
<b>Topband /attachment</b>	Silicone dot topband
<b>Key words</b>	Clima Comfort, Clima Fresh
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\* Made-to-measure products and variants please see our product catalogue.

# Measuring chart

mediven® 550 leg

Standard size	IV	V	VI	VII	VIII	
<p>l=Lengths    c=Circumferences</p>	<b>cD</b>	29–33	30–34	31–35	32–36	33–37
	<b>cC</b>	30–34	31–35	32–36	33–37	34–38
	<b>cB<sup>1</sup></b>	24–26	25–27	26–28	27–29	28–30
	<b>cB</b>	19–20	20–21	21–22	22–23	23–24
	<b>cY</b>	27–29	28–30	29–32	29–33	30–34
	<b>cA</b>	18–20	19–21	20–22	21–23	22–24
<b>Circumferences in cm</b>	<b>X</b>	<b>XI</b>	<b>XII</b>	<b>XIII</b>	<b>XIV</b>	
	<b>cD</b>	35–39	36–40	37–41	38–42	39–43
	<b>cC</b>	36–40	37–41	38–42	39–43	40–44
	<b>cB<sup>1</sup></b>	29–33	30–34	31–35	32–36	33–37
	<b>cB</b>	25–26	26–27	27–28	28–29	29–30
	<b>cY</b>	32–36	33–37	34–38	35–39	36–40
<b>cA</b>	24–26	25–27	26–28	27–29	28–30	
<b>Lengths</b>	<b>lD standard</b>	39–41 cm				
	<b>lD short</b>	37–39 cm				

# mediven mondi®


flat-knitted

<b>Product</b>	<b>mediven mondi*</b> Soft. Fine. Looks great.
<b>CCL</b>	2
<b>Lengths</b>	2 lengths for AD / calf and AG / thigh
<b>Standard colour</b>	Caramel
<b>Sizes</b>	I – VII
<b>Standard /made-to-measure</b>	Standard and made-to-measure
<b>Foot part</b>	Open toe straight ending
<b>Topband /attachment</b>	Silicone dot topband
<b>Key words</b>	Clima Comfort, Clima Fresh
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\* Made-to-measure products and variants please see our product catalogue.

# Measuring chart

mediven mondi®

Standard size	I	II	III	IV	V	VI	VII	
<b>ℓ=Lengths c=Circumferences</b> 	<b>cG with topband</b>	44–47	46–50	50–54	53–58	56–61	60–65	64–70
	<b>cG</b>	44–53	46–57	50–61	53–65	56–70	60–74	64–80
	<b>cF</b>	40–47	42–51	45–55	48–58	51–62	54–67	58–72
	<b>cE</b>	30–37	33–40	36–42	38–44	40–46	42–48	44–50
	<b>cD</b>	28–33	30–35	32–37	34–39	36–41	38–43	40–45
	<b>cC</b>	29–34	31–36	33–38	35–40	37–42	39–44	41–46
	<b>cB<sub>1</sub></b>	24–26	26–28	28–31	30–33	32–35	34–37	36–39
	<b>cB</b>	18–20	20–22	22–24	24–26	26–28	28–30	30–32
	<b>cY</b>	28–30	29–32	31–34	33–36	34–38	36–40	37–41
	<b>cA</b>	18–21	20–23	22–25	24–27	26–29	28–31	30–33
<b>Lengths</b>	<b>ℓG standard</b>	72–83 cm						
	<b>ℓG short</b>	62–71 cm						
	<b>ℓD standard</b>	38–43 cm						
	<b>ℓD short</b>	33–37 cm						

mediven esprit®

flat-knitted

<b>Product</b>	<b>mediven esprit*</b> Soft. Fine. Looks great.
<b>CCL</b>	1, 2
<b>Lengths</b>	Standard length: AG, CG Long length: AG, CG AC <sup>1</sup> hand piece
<b>Standard colour</b>	Caramel
<b>Sizes</b>	I – V
<b>Standard / made-to-measure</b>	Standard and made-to-measure
<b>Topband / attachment</b>	Narrow silicone dot topband
<b>Key words</b>	Clima Comfort, Clima Fresh
<b>Care instructions</b>	
<b>Material composition</b>	Elastane, Polyamide

\* Made-to-measure products and variants please see our product catalogue.

## Variants

### Styles


CG / SN

CG / SN long

AG / SN

AG / SN long

AC<sub>1</sub>



c = circumferences

## Standard sizes

Circumferences in cm	Sizes	I	II	III	IV	V
	cA		17.0 – 19.0	18.0 – 20.0	19.0 – 22.0	21.0 – 24.0
cC		15.0 – 17.0	17.0 – 19.0	19.0 – 21.0	21.0 – 23.0	23.0 – 25.0
cE		25.0 – 27.5	27.0 – 29.5	29.0 – 32.0	31.0 – 34.0	33.0 – 37.0
cG		28.0 – 32.0	30.0 – 34.0	32.0 – 37.0	35.0 – 40.0	38.0 – 44.0

### Lengths

CG standard 40 – 44 cm

CG long 44 – 48 cm

## Measuring chart

mediven esprit®

Circumference / Hand part (cm)					
Sizes	I	II	III	IV	V
a	17 – 20	19 – 22	21 – 24	23 – 26	25 – 28
b	18 – 21	20 – 23	22 – 25	24 – 27	23 – 25
c	15 – 17	17 – 19	19 – 21	21 – 23	23 – 25
c <sub>1</sub>	15.5 – 18	17.5 – 20	19.5 – 22	21.5 – 24	23 – 26

Length / Hand part (cm), all sizes	
AB	3.5 – 4.5
AC	10.0 – 12.0
AC <sub>1</sub>	15.0 – 16.0

## Measuring chart

mediven esprit®

Circumference / Fingers (cm)						
Variant AC1		I	II	III	IV	V
Thumb	top	4.5	5.5	6.5	7.5	8.5
	bottom	6.5	7.5	8.5	9.5	10.5
Index finger	top	4.5	5.5	6.5	7.5	8.5
	bottom	5.5	6.5	7.5	8.5	9.5
Middle finger	top	4.5	5.5	6.5	7.5	8.5
	bottom	5.5	6.5	7.5	8.5	9.5
Ring finger	top	4.0	5.0	6.0	7.0	8.0
	bottom	5.0	6.0	7.0	8.0	9.0
Little finger	top	3.5	4.5	5.5	6.5	7.5
	bottom	4.5	5.5	6.5	7.5	8.5

Length / Fingers (cm)	
Thumb	3.5 – 4.5
Index finger	4.5 – 5.5
Middle finger	5.5 – 6.5
Ring finger	4.5 – 5.5
Little finger	3.5 – 4.5

# mediven® standard colours

standard

Round-knitted Models	mediven elegance	mediven for men	mediven comfort	mediven plus	mediven active	mediven forte	mediven angio
Black	●	●	●	●	●	●	●
Anthracite	●	●	●	●	●	●	
Navy	●	●	●	●	●	●	
Cashmere	●		●	●		●	
Caramel	●		●	●		●	●
Rosé	●		○	●		○	
Bronze	●		○	●		○	
Beige	●		○	●		○	
medi Magenta	○			○			
White	○			○	●		
Brown		●			●		
Flat-knitted	mediven 550 leg	mediven mondi	mediven esprit	mediven harmony	<b>Note:</b> Made-to-measure products and variants are also available in trend colours, with Design- or Fashion-Elements and can also be embellished with crystal. For a complete overview please see our product catalogue.		
Sand	○	○	○	●			
Caramel	●	●	●	●			
Black	○	○	○	●			
Cashmere	○	○	○				
Anthracite	○	○	○				
Navy	○	○	○				

○ Delivery period up to 10 working days.